

Ascot

Rediscover the joy of mobility with the Ascot Mobility Scooter, your perfect companion for everyday adventures! Designed for effortless maneuverability, the Ascot glides smoothly through pavements, shopping centres, and local parks. Its lightweight structure, the heaviest part weighing just 30 lb (14 kg), makes it a breeze to lift and transport in your car, reigniting your freedom.

Comfort is at the heart of the Ascot's design. With a luxurious padded seat, you can enjoy lengthy outings without discomfort. The innovative tool-free disassembly allows you to break the scooter into five easy-to-handle parts, making loading and unloading a quick and simple task.

Safety and visibility are paramount, especially during those darker evenings. The Ascot is equipped with a bright headlight and luminous side lights, ensuring you are seen while you navigate your surroundings.

SPECIFICATION

Class	1	Total Weight	110lbs / 50kg
Length	40.1" / 102cm	Speed	4mph / 6kmh
Width	22" / 56cm	Battery	2x 12v 12Ah
Front	6.7" Solid	Capacity	Up to 21st 6lbs / 134.6kg
Rear	7.5" Solid	Max Slope	10.5% / 6°
Turning Radius	45.5" / 116cm	Lights	Front & Side
Ground Clearance	2.3" / 6cm	Suspension	No
Seat	Navigator		

COLOURS Red

Our policy is one of continual improvement.

We reserve the right to change specification without notice.



FreeRider®