

# Luggie

## Travel Chair

Sharing many benefits and design features of our original Luggie, our TravelChair is a smart, compact wheelchair which is easily foldable in just three steps. Equipped with joy-stick controls to offer an extra level of comfort as well as a soft seat cushion & solid tyres.

## Smart fold, Easy go

*joy-stick controls  
smart compartment  
extra comfortable chair*



## Specification

**Max Load:** 145kg / 320lbs  
**Total Weight:** 43kg / 90lbs  
**Battery Weight:** 2kg / 4.4lbs  
**Range:** up to 12km / 7.5miles\*  
**Max Speed:** 6.4km / 4mph  
**Max Incline:** 6 degrees  
**Turning Radius:** 77.9cm / 30.7inches  
**Ground Clearance:** 8 cm / 3.1 inches  
**Folded Dimensions:**  
45.7 x 60.6 x 95cm / 18 x 23.86 x 37.44inches (LxWxH)  
**Assembled Dimensions:**  
95 x 60.6 x 84cm / 37.445 x 23.86 x 33.11inches (LxWxH)

\* range varies depending on user weight & terrain. Based on 10.5ah lithium battery pack.  
optional upgrade to 16.5ah lithium battery for up to 12 miles range.

# FreeRider®