

Luggie

Elite

The Luggie Elite's patented design is the world's best selling foldable mobility scooter. Having changed the industry with its award winning, unique folding design which is often imitated, but never bettered! Rediscover your independence with a clear LED battery display, height adjustable seat & tiller column ensuring you are always riding in comfort & style.

Smart fold, Easy go

*extendable handlebars
clear LED battery display
adjustable armrests*



Specification

Max Load: 145kg / 320lbs
Total Scooter Weight: 25.8kg / 56.76 lbs
Battery Weight: 2kg / 4.4lbs
Range: up to 16km / 10miles*
Max Speed: 6.4km / 4mph
Max Incline: 6 degrees
Turning Radius: 103.8 cm / 40.9 inches
Ground Clearance:
Rear (at Motor) 3.7cm / 1.4 inches
Folded Dimensions:
39.5 x 45.5 x 66 cm / 15.5 x 17.91 x 26 inches (LxWxH)
Assembled Dimensions:
98.2 x 45.5 x 86 cm / 38.65 x 17.91 x 34 inches (LxWxH)

range varies depending on user weight & terrain. Based on 10.5ah lithium battery pack.

* optional upgrade to 16.5ah lithium battery for up to 12 miles range.

FreeRider®