

Luggie

Super Plus (3W)

Sharing many benefits and design features of our standard Super, yet exuding luxury and improved technology providing an enhanced user experience. Fitted with front wheel suspension and speed resistor technology automatically slowing as you turn corners, giving you peace of mind that you are riding safely.

Smart fold, Easy go

*extendable handlebars
additional under seat storage bag
wider seat base for added comfort*

Specification

Max Load: 163.6kg / 360lbs

Total Scooter Weight: 30.3kg / 66.7lbs

Battery Weight: 2kg / 4.4lbs

Range: up to 16km / 10miles*

Max Speed: 6km / 4mph

Max Incline: 6 degrees

Turning Radius: 95.25cm / 37.5inches

Ground Clearance:

Rear (at Motor) 4cm / 1.5inches

Folded Dimensions:

47.5 x 56.5 x 68cm / 18.7 x 22.4 x 26.7inches (LxWxH)

Assembled Dimensions:

99.5 x 56.5 x 83-88cm / 39.17 x 22.4 x 32.6-34.6inches (LxWxH)



range varies depending on user weight & terrain. Based on 10.5ah lithium battery pack.

* optional upgrade to 16.5ah lithium battery for up to 12 miles range.

FreeRider®