

# Luggie

## Elite Deluxe (4W)

Sharing many benefits and design features of our original Elite, yet exuding luxury and improved technology providing you with an enhanced user experience. A wider front wheelbase with added suspension, fitted with larger rear, solid tyres improving stability and comfort over those bumpier surfaces.

## Smart fold, Easy go

*height & width adjustable handlebars  
extra wide front wheel base  
front suspension*

### Specification

**Max Load:** 145kg / 320lbs

**Total Scooter Weight:** 28.8kg / 63.4lbs

**Battery Weight:** 2kg / 4.4lbs

**Range:** up to 16km / 10miles\*

**Max Speed:** 6.4km / 4mph

**Max Incline:** 6 degrees

**Turning Radius:** 92.5 cm / 36.4inches

**Ground Clearance:**

*Rear (at Motor) 4.5 cm / 1.77 inches*

**Folded Dimensions:**

99.5 x 51.5 x 39cm / 39.17 x 20.3 x 15.35inches (LxWxH)

**Assembled Dimensions:**

99.5 x 51.5 x 83-88cm / 39.17 x 20.3 x 32.6-34.6inches (LxWxH)



range varies depending on user weight & terrain. Based on 10.5ah lithium battery pack.

\* optional upgrade to 16.5ah lithium battery for up to 16 miles range.

# FreeRider®