

Luggie

Elite Plus (3W)

Sharing many benefits and design features of our original Elite, yet exuding luxury and improved technology providing an enhanced user experience. Fitted with speed resistor technology, suspension on the front wheel and larger rear solid tyres improving comfort and stability over those bumpier surfaces.

Smart fold, Easy go

*extendable handlebars
solid, puncture proof tyres
front end suspension*



Specification

Max Load: 145kg / 320lbs

Total Scooter Weight: 26.2kg / 57.6lbs

Battery Weight: 2kg / 4.4lbs

Range: up to 16km / 10miles*

Max Speed: 6.4km / 4mph

Max Incline: 6 degrees

Turning Radius: 95.25 cm / 37.5 inches

Ground Clearance:

Rear (at Motor) 4.5 cm / 1.77 inches

Folded Dimensions:

42.5 x 50 x 66 cm / 16.7 x 19.6 x 25.5 inches (LxWxH)

Assembled Dimensions:

99.5 x 50 x 83-88cm / 39.17 x 19.6 x 32.6-34.6inches (LxWxH)

* range varies depending on user weight & terrain. Based on 10.5ah lithium battery pack.

† optional upgrade to 16.5ah lithium battery for up to 16 miles range.

FreeRider®