

Luggie

Super Deluxe (4W)

Sharing many benefits and design features of our standard Super, yet exuding luxury and improved technology providing an enhanced user experience. Fitted with front wheel suspension and speed resistor technology, automatically slowing as you turn corners, giving you peace of mind that you are riding safely.

Smart fold, Easy go

*extendable handlebars
additional under seat storage bag
greater stability*

Specification

Max Load: 163.6 kg / 360 lbs

Total Scooter Weight: 32.1 kg / 70.6 lbs

Battery Weight: 2 kg / 4.4 lbs

Range: up to 16 km / 10 miles*

Max Speed: 6 km / 4 mph

Max Incline: 6 degrees

Turning Radius: 93.5 cm / 36.8 inches

Ground Clearance:

Rear (at Motor) 4 cm / 1.5 inches

Folded Dimensions:

99.5 x 56.5 x 41 cm / 39.17 x 22.4 x 16.1 inches (LxWxH)

Assembled Dimensions:

99.5 x 56.5 x 86-92.5 cm / 39.17 x 22.4 x 33-37.1 inches (LxWxH)



range varies depending on user weight & terrain. Based on 10.5ah lithium battery pack.

* optional upgrade to 16.5ah lithium battery for up to 12 miles range.

FreeRider®